

Docket No. MWH-0029US

PATENT

## IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

In re Application of: Stephen B. Liggett  
Application No.: 09/856,803  
Filed: May 25, 2001 (35 U.S.C. § 371 of PCT/US99/27963, filed November 24, 1999, which claims benefit of U.S. Appl. No. 60/109,886, filed November 25, 1998)  
Confirmation No.: 3706  
Group No.: 1634  
Examiner: Myers, C.  
For: **POLYMORPHISMS IN THE 5' LEADER CISTRON OF THE  $\beta_2$ -ADRENERGIC RECEPTOR**

Commissioner for Patents  
Washington, D.C. 20231

## Certificate of Facsimile Transmission

I hereby certify under 37 C.F.R. § 1.8 that this correspondence is being transmitted by facsimile to the United States Patent and Trademark Office, Commissioner for Patents, TC 1600, at (703) 872-9306, on May 7<sup>th</sup>, 2003.

  
Matthew M. Catlett

**DECLARATION OF STEPHEN B. LIGGETT, M.D., UNDER 37 C.F.R. § 1.131**

This Declaration Of Stephen B. Liggett, M.D., Under 37 C.F.R. § 1.131 is being submitted as part of Applicant's Response To Office Action Under 37 C.F.R. § 1.11 regarding the office action dated January 7, 2003 that was received in the captioned application.

Being warned that willful false statements and the like are punishable by fine or imprisonment, or both, under 18 U.S.C. § 1001, and that such willful false statements and the like may jeopardize the validity of the instant application or patent resulting therefrom, I hereby declare that:

1) I am the original, sole, and first inventor of the subject matter that is claimed in

pending claims 1-8 and 11 of the captioned application, namely (a) a method for removing the

Docket No. MWH-0029US

PATENT

$\beta_2$ -adrenergic receptor ( $\beta_2$ AR) gene of an individual comprising determining the identity of the nucleotide pair at the 5' leader cistron (5'LC) polymorphic site (PS), which, as is demonstrated throughout the specification of the captioned application, is located 47 bases upstream of the  $\beta_2$ AR coding region, which begins at nucleotide position 1588 of SEQ ID NO:1 (thus, the 5'LC PS is located at nucleotide position 1541 of SEQ ID NO:1) in the two copies of the  $\beta_2$ AR gene present in the individual; and (b) a method for genotyping the  $\beta_2$ AR gene of an individual comprising determining the identity of the nucleotide pair at the 5'LC PS and at one or more additional PSs in the  $\beta_2$ AR gene in the two copies of the  $\beta_2$ AR gene present in the individual.


2) Further to an effort, dating back to as early as January of 1996 (see attached copies of PCR protocols), to discover polymorphisms in the region upstream of the  $\beta_2$ AR gene, I directed the performance of an experiment designed to elucidate the existence, if any, of such polymorphisms. Utilizing PCR techniques to analyze genomic DNA in this region from human volunteers, I discovered, in the "sense" strand, the existence of a thymine residue 47 bases upstream of the  $\beta_2$ AR coding region, as well as the existence of an adenine residue 47 bases upstream of the  $\beta_2$ AR coding region in the "antisense" strand. Copies of chromatograms generated by the automated sequencer used to sequence the PCR products demonstrating this discovery are attached (chromatogram #096-1369 demonstrates a thymine residue in the "sense" strand at the nucleotide position that is located 47 bases upstream of the  $\beta_2$ AR coding region; chromatogram #096-1364 demonstrates an adenine residue in the "antisense" strand at the nucleotide position that is located 47 bases upstream of the  $\beta_2$ AR coding region; chromatogram #096-1367 demonstrates a thymine residue in the "sense" strand at the nucleotide position that is located 47 bases upstream of the  $\beta_2$ AR coding region; and chromatogram #096-1362 demonstrates an adenine residue in the "antisense" strand at the nucleotide position that is located 47 bases upstream of the  $\beta_2$ AR coding region). Although all previous reports indicated that the only known residue at the nucleotide position located 47 bases upstream of the  $\beta_2$ AR

Docket No. MWH-0029US

PATENT

coding region, in the "sense" strand, was a cytosine (and thus, in the "antisense" strand, a guanine), to confirm that I had indeed discovered a polymorphism at this position, I subsequently directed the performance of a similar experiment with the wild-type sequence, and discovered, in the "sense" strand, a cytosine, and in the "antisense" strand, a guanine. Copies of chromatograms generated by the automated sequencer used to sequence the PCR products demonstrating this discovery are attached (chromatogram #096-2859 demonstrates a cytosine residue in the "sense" strand at the nucleotide position that is located 47 bases upstream of the  $\beta_2$ AR coding region; and chromatogram #096-2860 demonstrates a guanine residue in the "antisense" strand at the nucleotide position that is located 47 bases upstream of the  $\beta_2$ AR coding region). My discovery of this polymorphism, and my subsequent confirmation of this discovery, occurred prior to the effective date of any of the following references: Timmermann *et al.*, *Kidney Intl.* 53:1455-60 (June 1998), Timmerman *et al.*, *J. Molecular Med.* 76:B30, Abst. P-109 (May 1998), Timmermann *et al.*, *Human Mutation* 11(4):343-4 (March 1998). With respect to the copies of the chromatograms, the nucleotide position that is 47 bases upstream of the  $\beta_2$ AR coding region is that denoted with a "^" symbol.

3) All statements made herein of my knowledge are true, and all statements made herein on information and belief are believed by me to be true.



Stephen B. Liggett, M.D.  
Director, Division of Pulmonary and  
Critical Care Medicine  
University of Cincinnati Medical Center








七


1-3231

C N C N T C T A G G G C C C A G C C A G G T A G C G C G G G A A G C A T A G T G G C C C C C C C C T C A G G G M A S C A G T T A G A C C C G G C C C C C C G A C C A C C O T C A G A N A N A G A A G G C G N A G G A N G A G T A A


NO. 730




120  
 130  
 140  
 150  
 160  
 170  
 180  
 190  
 200  
 210  
 220



CAGAAAGGACTGCGAAGCGGCTTCTGGAGACACGGGTGGAACCTGGCAGCACCCGCGAAGCCCCATAGCAACCACAAGTGAGTGTGGAGAGAAGANNTGCCAACCAAC



CCGACAGCAGCGGGCTGAATGAAAGCTTCGAGAGCGTCCGTCTGCAGGCCGCAAAACCCTGCCCATGGATACCTGTTGCTGTAAGCCCCAATGCGCT

[illegible]

5. MINN CO

VAY. 7. 2003 6:27PM

GENAISANCE PHARM.








Spacing: 10.31

3

UACAAAGGA GGGG GGCG TGGCGGCTTGGCCATGCCCCGCACAT CTGGAGGTGAGCGCACT GECT AGGGGAGGCTCAGCGGAGGAGAGCCACAGGCGGGGCT CTACGAGGCCGGA A

P. 132

MAY. 7. 2003 6:28PM



PCR (B2AR)

1/26/96

\* Made new set of primers for B2AR 5' flanking region. ~~both~~  
Product should include the CRE + short ORF lines  
were chosen = MacVector

blood samples are being obtained from pts in the UC cell line  
clinic by Melanie Meyers.

DNA samples prepared by big Donnelly using Dinitrogen + 10  
Dianc kits. Samples are numbered as received (A1, A2, etc.)

- Set up PCR rxn using new primers + made up 96ul master mix.

10ul buffer II

6ul 25mM MgCl<sub>2</sub>

0.8ul 25mM dNTP

0.5ul 100uM forward primer (B2AR-F1)

0.5ul 100uM reverse primer (B2AR-R1)

77.7ul dH<sub>2</sub>O

0.5ul Taq polymerase

- digest 24ul of master mix into 8 PCR tubes

- add 1ul template DNA

- overlay = 1 drop mineral oil

- perform PCR in thermocycler using the following gradient made on


94°C x 2 min

94°C x 30 sec / 64°C, or 62°C, 60°C, or 58°C x 30 sec / 72°C x 30 sec => 35 cycles

72°C x 7 min

- run 10ul of PCR rxn on TAE minigel

1/26/96 (cont)



1 - 100 bp ladder

6 - blank

2 - A1 (Amplify) 64°C

7 - A1 (Amplify) 64°C

3 - " 62°C

8 - " 62°C

4 - " 60°C

9 - " 60°C

5 - " 58°C

10 - " 58°C

No product seen = DNA isolated by Amplify kit. Nice band seen = DNA from Qiana kit but expected size should be 528 bp. The band present appears to be < 400 bp.


Repeat PCR run with following modifications:

- add 73.7 ul dH<sub>2</sub>O

- aliquot 23 ul of master mix into PCR tubes

- add 2 ul template DNA

- change cycling temp to Amplify template to 56°, 54°, 52°, 50°C; change cycling temp for Qiana template to 58°, 56°, 54°, 52°C.



1 - 100 bp ladder

6 - 100 bp ladder

2 - A1 (Amplify) 56°C

7 - A1 (Amplify) 56°C

3 - " 54°C

8 - " 56°C

4 - " 52°C

9 - " 54°C

5 - " 50°C

10 - " 52°C

Now have band in all samples. However, uppermost band present in both samples sets appears to be < 500 bp, still smaller than expected.

PCR (5 min)

Set up master mix for four (4) 25ul PCR rxns:

2 ul template DNA (DWM)  
 1.5 ul forward primer  
 1.5 ul reverse primer  
 20 ul 5X buffer (buffer A from Stratagene PCR optimization kit)  
 10 ul dUTP (2.5mM) (from optimization kit)  
 65 ul dH<sub>2</sub>O  
0.8 ul tag  
 100 ul total

- aliquot 25 ul of master mix into 4 PCR rxn tubes
- overlay  $\bar{c}$  : drop of mineral oil
- PCR cycle
 

98°C	2 min	} 30 cycles
98°C	30 sec	
56° 54° 52° ~ 50°	30 sec	
72°C	30 sec	
72°C	7 min	

- remove 10ul aliquot &amp; run on minigel

#1 - 100 bp ladder  
 #2 - 56°C  
 #3 - 54°C  
 #4 - ~~56~~ 50°C  
 #5 - 50°C

